

# Level 2 Protective Gown

Personal Protective Equipment

## Features

- Single-use Disposable
- Tear Resistant
- Breathable AP Coated
- Protects Against seepage
- Anti-Static
- Non-Surgical
- Dirt Resistant
- Made from Non-Woven Polyethylene Fiber
- FOB CA Pricing CNS 14798:2004,
- Certification CNS 12915:2010, CNS 10460:2007,  
CNS 14801:2004, CNS 8150:1981,  
CNS 12222:2009



## Scope of Use

This Level 2 Protective Gown provides bacterial and seepage-proofing, and has been shown to be effective during low-risk, general medical, and non-surgical use.

**Package:** Individually

**HCA Part Number:** OP-IG-L2-EME

## Storage

Level 2 Gowns are designed to protect users so its important to store them as follows:

- In original package
- In a clean air environment free from hazards
- Away from direct sunlight
- In a cool temperature controlled environment

## Instructions For Use

Put on Level 2 Protective gown and make sure flaps are closed properly. Use adhesive tapes to seal off opening around face and wrists to maximize protection.